




Agoralaan Abis • 3590 Diepenbeek • Belgium
Tel: +32(0)11/28.69.02 • www.bedal.be

BEDAL PRODUCT CATALOGUE

Version V1, June 2019.

bedalpatch

IV shower protection dressing. The Bedalpatch is a single use dressing for covering an intravenous catheter and its exit site during washing and showering. The Bedalpatch is available in 2 sizes: a large size for the chest position and a small size for the arm position.




Product name	Product code	Classification	Description	# per box	# box per case
BEDALPATCH®	BD001	Class I non sterile	Shower dressing large for chest position	100	1
BEDALPATCH® comfort	BD003	Class I non sterile	Shower dressing small for arm position, PICC	100	1

drainpatch

Catheter securement device. The Drainpatch is a sterile, single use patch for the stabilization of catheters. It can be applied for catheters with a French size between 6Fr and 16Fr.

The Drainpatch is available sterile and individually pouched, or as a bulk non sterile to include in procedure packs. For inclusion in procedure packs a validation of this setup is required.



Product name	Product code	Classification	Description	# per box	# box per case
DRAINPATCH®	DP01	Class I sterile	Catheter stabilization device for 6-16 French, sterile	25	4
DRAINPATCH® bulk non-sterile	DPNS01	Class I sterile	Catheter stabilization device for 6-16 French, bulk non-sterile	25	4